

HCR RNA-FISH/IF for samples on slide

This protocol has not been validated for all mammalian cell types and should only be used as a template.
Technical support: support@moleculartechnologies.org

Sample preparation

Samples should be prepared in the same manner as for a traditional immunohistochemistry, up to the primary antibody application step. This may include permeabilization and antigen retrieval. Then proceed with the protocol described below.

HCR IF with/without HCR RNA-FISH

Unlabeled primary antibody probes and initiator-labeled secondary antibody probes for protein targets, split-initiator DNA probes for RNA targets, and simultaneous HCR signal amplification for all targets.

Protein detection stage

1. Block samples with 200 μL of antibody buffer on top of the sample. Incubate at room temperature for 1 h in a humidified chamber.
2. Prepare working concentration of primary antibodies in antibody buffer. Prepare 100 μL per chamber.
NOTE: follow manufacturer's guidelines for primary antibody working concentration.
3. Drain slide by blotting edges on a Kimwipe and wipe around the section with another Kimwipe.
4. Add primary antibody solution to each section and incubate overnight (>12 h) at 4 °C in a humidified chamber.
NOTE: Incubation may be optimized (e.g., 1–2 h at room temperature) depending on sample type and thickness.
5. Remove excess antibodies by immersing slide in 1× PBST at room temperature 3 × 5 min.
6. Prepare 1 $\mu\text{g}/\text{mL}$ working concentration of initiator-labeled secondary antibodies in antibody buffer. Prepare 100 μL per chamber.
NOTE: Concentration may be optimized depending on protein target and primary antibody.
7. Drain slide by blotting edges on a Kimwipe and wipe around the section with another Kimwipe.
8. Add secondary antibody solution to each section and incubate for 1 h at room temperature in a humidified chamber.
9. Remove excess antibodies by immersing slide in 1× PBST at room temperature 3 × 5 min.
10. Proceed to **RNA detection stage** for co-detection of protein and RNA. Otherwise, proceed to **Amplification stage**.

RNA detection stage

1. Drain slide by blotting edges on a Kimwipe and wipe around the section with another Kimwipe.
2. Post-fix sample with 200 μL of 4% formaldehyde.
CAUTION: use formaldehyde with extreme care as it is a hazardous material.
3. Incubate for 10 min at room temperature.
4. Immerse slides for 2×5 min in PBST.
5. Immerse slides for 5 min in $5 \times$ SSCT.
6. Pre-warm a humidified chamber to 37°C .
7. Drain slide by blotting edges on a Kimwipe and wipe around the section with another Kimwipe.
8. Pre-hybridize samples in 200 μL of probe hybridization buffer for 10 min at 37°C .
CAUTION: Probe hybridization buffer contains formamide, a hazardous material.
NOTE: pre-heat probe hybridization buffer to 37°C before use.
9. Prepare a 16 nM probe solution by adding 1.6 pmol of each probe mixture (e.g. 1.6 μL of 1 μM stock) to 100 μL of probe hybridization buffer at 37°C .
NOTE: This is the amount of probe set needed for each target on a single sample using 100 μL of incubation volume.
10. Remove the pre-hybridization solution and drain excess vuffer on slide by blotting edges on a Kimwipe.
11. Add 100 μL of the probe solution on top of the tissue sample.
12. Place a coverslip on the sample and incubate overnight (>12 h) in the 37°C humidified chamber.
13. Immerse slide in probe wash buffer at 37°C to float off coverslip
CAUTION: Probe wash buffer contains formamide, a hazardous material.
NOTE: pre-heat probe wash buffer to 37°C before use.
14. Remove excess probes by incubating slides at 37°C
 - (a) 75% of probe wash buffer / 25% $5 \times$ SSCT for 15 min
 - (b) 50% of probe wash buffer / 50% $5 \times$ SSCT for 15 min
 - (c) 25% of probe wash buffer / 75% $5 \times$ SSCT for 15 min
 - (d) 100% of probe wash buffer / 75% $5 \times$ SSCT for 15 min
15. Proceed to **Amplification stage**.

Amplification stage

1. Immerse slide in $5\times$ SSCT at room temperature for 5 min.
2. Drain slide by blotting edges on a Kimwipe and wipe around the section with another Kimwipe.
3. Pre-amplify samples in $200\ \mu\text{L}$ of amplification buffer in a humidified chamber for 30 min at room temperature.
NOTE: Equilibrate amplification buffer to room temperature before use.
4. Separately prepare 6 pmol of hairpin h1 and 6 pmol of hairpin h2 by snap cooling $2\ \mu\text{L}$ of $3\ \mu\text{M}$ stock (heat at $95\ ^\circ\text{C}$ for 90 seconds and cool to room temperature in a dark drawer for 30 min).
NOTE: HCR hairpins h1 and h2 are provided in hairpin storage buffer ready for snap cooling. h1 and h2 should be snap cooled in separate tubes. This is the amount of hairpins needed for each target in a single sample using $100\ \mu\text{L}$ of incubation volume.
5. Prepare a 60 nM hairpin solution by adding all snap-cooled h1 hairpins and snap-cooled h2 hairpins to $100\ \mu\text{L}$ of amplification buffer at room temperature per sample.
6. Remove the pre-amplification solution and drain excess buffer on slide by blotting edges on a Kimwipe.
7. Add $100\ \mu\text{L}$ of the hairpin solution on top of the tissue sample.
8. Incubate the slide overnight ($>12\ \text{h}$) in a dark humidified chamber at room temperature.
9. Remove excess hairpins by immersing slides in $5\times$ SSCT at room temperature for:
 - (a) $2\times 5\ \text{min}$
 - (b) $2\times 15\ \text{min}$
 - (c) $1\times 5\ \text{min}$
10. Drain slide by blotting edges on a Kimwipe and dry around the section with another Kimwipe.

Sample mounting for microscopy

1. Remove final wash and add $50\text{--}100\ \mu\text{L}$ of antifade mounting medium (e.g., Fluoromount-G with DAPI).
2. Place a coverslip on top for microscopy.
3. Slides can be stored at $4\ ^\circ\text{C}$ protected from light prior to imaging.

Buffers for HCR IF with/without HCR RNA-FISH

HCR probes (initiator-labeled antibody probes, split-initiator DNA probes), amplifiers, and buffers (antibody buffer, probe hybridization buffer, probe wash buffer, amplification buffer) are available from Molecular Technologies (www.moleculartechnologies.org). Probe hybridization buffer and probe wash buffer should be stored at -20 °C. Antibody buffer and amplification buffer should be stored at 4 °C. Make sure all solutions are well mixed before use.

1× PBST

1× phosphate-buffered saline (PBS)
0.1% Tween 20

For 40 mL of solution

4 mL of 10× PBS
400 μ L of 10% Tween 20
Fill up to 40 mL with ultrapure H₂O

5× SSCT

5× saline sodium citrate (SSC)
0.1% Tween 20

For 40 mL of solution

10 mL of 20× SSC
400 μ L of 10% Tween 20
Fill up to 40 mL with ultrapure H₂O

NOTE: avoid using calcium chloride and magnesium chloride in PBS as this leads to increased autofluorescence in the samples.

Reagents and supplies

ibidi μ -slide ibitreat (ibidi Cat. # 80826)
Poly-D-lysine hydrobromide (Sigma-Aldrich Cat. # P7280)
Molecular biology grade H₂O (Corning Cat. # 46-000-CV)
DPBS, no calcium, no magnesium (Life Technologies Cat. # 14190144)
Image-iT Fixative Solution 4% (Thermo Fisher Scientific Cat. # FB002)
10× PBS (Ambion Cat. # AM9624)
10% Tween 20 (Teknova Cat. # T0710)
20× saline sodium citrate (SSC) (Life Technologies Cat. # 15557-044)
DAPI Fluoromount-G (SouthernBiotech Cat. # 0100-20)

Citation Notes

Molecular Technologies (moleculartechnologies.org) is a non-profit academic resource within the Beckman Institute at Caltech. For citation, please select from the list below as appropriate for your application:

- **10-plex HCR spectral imaging**
HCR spectral imaging of any combination of 10 RNA and/or protein targets with 1-step quantitative HCR signal amplification for all targets simultaneously (Schulte *et al.*, 2024a).
Software: [HCR Imaging Python module including Dot Detection 2.0 and Unmix 1.0 packages](#).
- **HCR protein:protein complex imaging**
Multiplexed, quantitative, high-resolution imaging of protein:protein complexes with/without HCR IF and/or HCR RNA-FISH with simultaneous HCR signal amplification for all protein:protein, protein, and RNA targets (Schulte *et al.*, 2024b).
- **HCR RNA-FISH/IF**
A unified approach to multiplexed, quantitative, high-resolution RNA fluorescence in situ hybridization (RNA-FISH) and protein immunofluorescence (IF), with quantitative 1-step enzyme-free HCR signal amplification performed for all RNA and protein targets simultaneously (Schwarzkopf *et al.*, 2021).
- **HCR IF**
Multiplexed, quantitative, high-resolution protein immunofluorescence (IF) in highly autofluorescent samples (e.g., FFPE brain tissue sections) (Schwarzkopf *et al.*, 2021).
- **HCR RNA-FISH (v3.0)**
Automatic background suppression for dramatically enhanced performance (signal-to-background and quantitative precision) and ease-of-use (no probe set optimization for new targets) (Choi *et al.*, 2018). v3.0 supports three quantitative modes:
 - **qHCR imaging:** analog mRNA relative quantitation with subcellular resolution in an anatomical context;
 - **qHCR flow cytometry:** analog mRNA relative quantitation for high-throughput single-cell analysis;
 - **dHCR imaging:** digital mRNA absolute quantitation with single-molecule resolution in an anatomical context.
 Protocols for v3.0 in diverse organisms are adapted from the zoo paper.
Software: [Dot Analysis 1.0 package](#).
- **Quantitative HCR (qHCR) imaging with subcellular resolution**
Analog mRNA relative quantitation with subcellular resolution in the anatomical context of whole-mount vertebrate embryos. The read-out/read-in analysis framework enables bidirectional quantitative discovery in an anatomical context (Trivedi *et al.*, 2018).
Software: [Read-out/Read-in 1.0 package](#).
- **Zoo paper**
Protocols for multiplexed mRNA imaging in diverse sample types (Choi *et al.*, 2016):

bacteria in suspension	whole-mount zebrafish embryos and larvae
whole-mount fruit fly embryos	whole-mount chicken embryos
whole-mount worm larvae	whole-mount mouse embryos
whole-mount sea urchin embryos	FFPE human tissue sections
- **Digital HCR (dHCR) imaging with single-molecule resolution**
Digital mRNA absolute quantitation with single-molecule resolution in thick autofluorescent samples (e.g., 0.5 mm adult mouse brain sections) (Shah *et al.*, 2016).

- **Multiplexed quantitative HCR (qHCR) northern blots**
Simultaneous quantification of RNA target size and abundance with signal amplification for up to 5 target RNAs (Schwarzkopf & Pierce, 2016).
- **HCR RNA-FISH (v2.0)**
2nd generation HCR RNA-FISH technology (v2.0) using DNA probes and DNA HCR amplifiers: 10× increase in signal, 10× reduction in cost, dramatic increase in reagent durability (Choi *et al.*, 2014).
- **HCR RNA-FISH (v1.0)**
1st generation HCR RNA-FISH technology (v1.0) using RNA probes and RNA HCR amplifiers: multiplexed mRNA imaging in whole-mount vertebrate embryos with simultaneous signal amplification for up to 5 target mRNAs (Choi *et al.*, 2010).
- **HCR mechanism**
The hybridization chain reaction (HCR) mechanism enables isothermal enzyme-free molecular signal amplification (Dirks & Pierce, 2004).

We are happy to provide advice and technical support: support@moleculartechnologies.org.

– The Molecular Technologies Team

References

- Choi, H. M. T., Chang, J. Y., Trinh, L. A., Padilla, J. E., Fraser, S. E., & Pierce, N. A. (2010). Programmable in situ amplification for multiplexed imaging of mRNA expression. *Nat Biotechnol*, **28**(11), 1208–12. ([pdf](#), [supp info](#), [supp movie 1](#), [supp movie 2](#))
- Choi, H. M. T., Beck, V. A., & Pierce, N. A. (2014). Next-generation in situ hybridization chain reaction: higher gain, lower cost, greater durability. *ACS Nano*, **8**(5), 4284–4294. ([pdf](#), [supp info](#), [supp movie](#))
- Choi, H. M. T., Calvert, C. R., Husain, N., Huss, D., Barsi, J. C., Deverman, B. E., Hunter, R. C., Kato, M., Lee, S. M., Abelin, A. C. T., Rosenthal, A. Z., Akbari, O. S., Li, Y., Hay, B. A., Sternberg, P. W., Patterson, P. H., Davidson, E. H., Mazmanian, S. K., Prober, D. A., van de Rijn, M., Leadbetter, J. R., Newman, D. K., Readhead, C., Bronner, M. E., Wold, B., Lansford, R., Sauka-Spengler, T., Fraser, S. E., & Pierce, N. A. (2016). Mapping a multiplexed zoo of mRNA expression. *Development*, **143**, 3632–3637. ([pdf](#), [supp info](#), [fruit fly movie](#), [sea urchin movie](#), [zebrafish movie](#), [chicken movie](#), [mouse movie](#))
- Choi, H. M. T., Schwarzkopf, M., Fornace, M. E., Acharya, A., Artavanis, G., Stegmaier, J., Cunha, A., & Pierce, N. A. (2018). Third-generation in situ hybridization chain reaction: multiplexed, quantitative, sensitive, versatile, robust. *Development*, **145**, dev165753. ([pdf](#), [supp info](#), [Dot Analysis 1.0 package](#))
- Dirks, R. M., & Pierce, N. A. (2004). Triggered amplification by hybridization chain reaction. *Proc Natl Acad Sci USA*, **101**(43), 15275–15278. ([pdf](#))
- Schulte, S.J., Fornace, M.E., Hall, J.K., Shin, G.J., & Pierce, N.A. (2024a). HCR spectral imaging: 10-plex, quantitative, high-resolution RNA and protein imaging in highly autofluorescent samples. *Development*, **151**(4), dev202307. ([pdf](#), [supp info](#), [HCR Imaging Python Module containing Dot Detection 2.0 and Unmix 1.0 packages](#))
- Schulte, S.J., Shin, B., Rothenberg, E.V., & Pierce, N.A. (2024b). Multiplex, quantitative, high-resolution imaging of protein:protein complexes via hybridization chain reaction. *ACS Chem Biol*, **19**(2), 280–288. ([pdf](#), [supp info](#))
- Schwarzkopf, M., & Pierce, N. A. (2016). Multiplexed miRNA northern blots via hybridization chain reaction. *Nucleic Acids Res*, **44**(15), e129. ([pdf](#), [supp info](#))
- Schwarzkopf, M., Liu, M. C., Schulte, S. J., Ives, R., Husain, N., Choi, H. M. T., & Pierce, N. A. (2021). Hybridization chain reaction enables a unified approach to multiplexed, quantitative, high-resolution immunohistochemistry and in situ hybridization. *Development*, **148**(22), dev199847. ([pdf](#), [supp info](#))
- Shah, S., Lubeck, E., Schwarzkopf, M., He, T.-F., Greenbaum, A., Sohn, C. H., Lignell, A., Choi, H. M. T., Gradi-naru, V., Pierce, N. A., & Cai, L. (2016). Single-molecule RNA detection at depth via hybridization chain reaction and tissue hydrogel embedding and clearing. *Development*, **143**, 2862–2867. ([pdf](#), [supp info](#), [supp movie 1](#), [supp movie 2](#), [supp movie 3](#))
- Trivedi, V., Choi, H. M. T., Fraser, S. E., & Pierce, N. A. (2018). Multidimensional quantitative analysis of mRNA expression within intact vertebrate embryos. *Development*, **145**, dev156869. ([pdf](#), [supp info](#), [Read-out/Read-in 1.0 package](#))